SHOCK  
DOCTOR

GEL MAX POWER

6900 - 6901 - 6954

EN	SPORTS MOUTHGUARD
FR	PROTEGE-DENTS DÉSPORT
DE	SPORTMUND SCHUTZ
IT	PARADENTOTICO SPORTIVO
ES	PROTECTOR BUCAL DEPORTIVO
PT	PROTECTOR BUCAL ESPORTIVO
NL	SPORTMOUTDENGUARD
BG	СТИХОТВОРЕНІЯ ПРОДАВЦІЯ ЗАСЛУГ
HR	SPORTSKA USTAVITNICKA
CS	SPORTOVNÍ ÚSTNÍ CHRÁNIČ

DA	SPORTMUNDVAGT
ET	SPORDI SUUKTASE
FI	URHEILUNSUUSKOJUS
EL	AEGAHITOKS PROTEKTÖÖRETZOMATO
HU	PROTEKTÖÖRGÖVÖL
LV	SPORTA MŪTĀS SĀRS
LT	SPORTINIS BURNOS APSAUGAS
PL	OCHRONAJAZIUS LUSTPORTOVY
RO	PROTECTIE BUCAL SPORTIVĂ
SK	SPORTOVNÝ ČERNÝ ČRHNÁČ
SL	SPORTNIČKI
SV	SPORTMUNSKYDD
NO	IDRETTSMUNSKJØLD

©2004 United Sports Brands  
11880 Valley Center, CA 92078  
1-800-233-6956  
For Info: [www.shockdoctor.com/](http://www.shockdoctor.com/)

Imported / Distributed by  
United Sports Europe  
Janssen-Pharmaceuticals 4  
Post Box 4  
2440 Gelgau, Belgium

UNITED SPORTS  
United Sports  
Sportsbrands  
Sportswear  
Westbahnhofstrasse 1  
CH-4500 Solothurn, Switzerland  
+32 (0) 14 85 47 85

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QNET LTD  
Livingstone House,  
303 Harrow Road  
Wembley, Middlesex HA9 6BD  
United Kingdom  
+44-2045250603

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DA BRUGERINFORMATION

**⚠️WARNING**

This mouthguard is engineered to mitigate dental injuries during specific sports and recreational activities like rugby, hockey, boxing, martial arts, and more. However, it is important to note that this mouthguard is not designed to prevent severe impact or prevent critical head injuries and falls. Incorrect molding or usage can compromise protection. This device is specifically designed for individual use and is not intended for multiple users. A mouthguard or use application incorrect may compromise the protection. Get your personal dentist or orthodontist to fit the mouthguard for you. It is also used for those conditions where the bruising inspection is required on any signs of wear or damage, and replace immediately if any extreme conditions, or using it beyond its intended scope significantly reduces its effectiveness and safety. Sharing or chewing on the mouthguard can also damage the mouthguard.

**PLEASE NOTE:** This mouthguard's fitting process requires boiling water, ensure to handle with care, and do not boil water for more than 10 minutes under 18, adult supervision is necessary to prevent injury.

**IMPORTANT WARNING:** This mouthguard should not be used for dental braces or other dental devices, unless approved by the person fitting the brace. You must consult a dentist or orthodontist who is fitting a brace or dental device to ensure compatibility and safety of this or any other mouthguard. It is essential to ensure compatibility and safety under such conditions. Do not allow anyone other than the user to handle or store the mouthguard. If anyone else handles or stores the mouthguard, then the user must be advised to replace it as soon as possible.

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